

Cerkamec

Made in Poland

MEDICAL COMPANY

imted Edition

New golden tip:












Join the world's endodontic elite with our a limited edition of BTR PEN














Dental tool for removing broken canal instruments.



Seize the impossible!

single-tool system for every type of root canal



precise grip of the broken instrument



thinnest working tip on the market - only Ø 0,2 mm



the body made of surgical steel




life-long warranty

Only one on the market!

The tool combines new technology, fine quality and effective work. BTR Pen makes the removal of broken dental instruments a predictable, universal and repeatable procedure.

The usage of the ultrathin and highly-elastic working tip, togehter with ultimate-endurance nitinol loop, makes it easy to get to the broken canal instrument. The tool allows certain grip and safe removal of the broken canal instrument regardless of it's possition in the canal.







Included:

Spring:

- full control of loop's tension and size
- optimalized elasticity parameter
- ya- guarantee of light and flawless moves
 - freedom of fitting the spring toughness

Changeable tip with nitinol loop:

- the combination of a thin and highly elastic needle of Ø 0,2 mm with nitinol loop
- easiness in insertion into narrow and bent canals
- primary shape memory
- maximum endurance to stretching
- biocompatibility, corrosion resistance



- ergonomic and comfortable handle equipped with grooves facilitating a secure hand grip
- made from high quality surgical stainless steel
- universal size
- designed for multiple use
- safe for autoclave sterilisation



- quick and easy formation of the loop into desired size
- made from stainless steel
- designed for multiple use



Impossible changes into possible

Broken Remover

New tips with new sizes now available

l mm

(Available soon)

Thinnest tip on the market!



NEW!

Ideal kit for training

the removal of broken endodontic files from root canals

The set includes:

btf (

Dental training

DIGITAL MICROSCOPE

ADSM201

- MINI HDMI input
- MINI USB port
- MICRO SD card input
- hight regulation
- LED lightning

MICROSCOPE'S INSTRUCTION FOR USE

Dental training penal peral pe

BTR Pen WORK METHODOLOGY

work with BTR Pen step by step

ANTI-SLIP MAT

· stability during the work with the kit

▼ BTR DENTAL TRAINING KIT

TWO TRAINING TEETH

with pieces of broken endodontic files

HARD 3 deeply placed, blocked pieces of files

(requiring processing before removal)

EASY 3 loosely placed pieces of files

SILICONE TRAINING MASS

- · transparent and backlit
- · 8 pieces of broken endodontic files
- · a space to place the training tooth

SILICONE MASS STAND

allows a 3D rotation

Also available:

Trainin h

The set includes:

SILICONE TRAINING MASS

with pieces of broken endodontic files

SILICONE MASS STAND

· allows a 3D rotation

ANTI-SLIP MAT

4

BTR DENTAL TRAINING GUM



easy

TWO TRAINING TEETH

with pieces of broken endodontic files (easy/hard)

10 ENDODONTICS





- Available packages: bottle of 200 g bottle of 400 g

- use with: ENDO-PREP CREAM









Available packages: - bottle of 200 g - bottle of 400 g

use with: ENDO-PREP CREAM







▲ CHLORAXID GEL EXTRA

Preparation for root canals rinsing

Sodium hypochlorite 2% or 5,25% in gel + surfactants









14 ENDODONTICS





Available packages: - bottle of 200 g - bottle of 400 g

use with: HYDROCAL



▲ GLUCO-CHEX 2%

Chlorhexidine digluconate 2% in liquid

Available packages: - bottle of 200 g - bottle of 400 g



CITRIC KWAS CYTRYNOWY 40% CYTRYNOWY 4

▲ CITRIC ACID 40% LIQUID for root canals rinsing





▲ ISOPROPYL ALCOHOL







Available for products:

CHLORAXID 2% EXTRA,

ENDO-SOLUTION PREMIUM. CITRIC ACID 40% 20 x 20 x

Available for products: CHLORAXID 2% EXTRA

CHLORAXID 2%,

CHLORAXID 5,25%,

GLUCO-CHEX 2%,

FNDO-SOLUTION

CHLORAXID 5,25% EXTRA,









The best in the world!



BEST

short setting

plasticity



Biomaterial for filling and rebuilding of root canals/ Biomaterial do wypełniania i odbudowy kanałów korzeniowych





Available packages:

- MINI 3 x 0,14 g BIO MTA+ powder, 1 ml of BIO MTA+ liquid
- STANDARD 6 x 0,14 g BIO MTA+ powder, 1 ml of BIO MTA+
- MAXI 10 x 0,14 g BIO MTA+ powder, 2 x 1 ml of BIO MTA+ liquid
- PRO 10 x 0,14 g BIO MTA+ powder, 2 x 1 ml of BIO MTA+ liquid, MATRIX MTA+ $\,$

biocompatibility

strength

BIO MTA+ A Biomaterial for filling and rebuilding of root canals



Available packages:

- MINI 3 x 0,14 g MTA+ powder, 1 ml liquid
- STANDARD 6 x 0,14 g MTA+ powder, 1 ml liquid
- MAXI 10 x 0,14 g MTA+ powder, 1 ml liquid
- PRO 10 x 0,14 g MTA+ powder, 1 ml liquid, MATRIX MTA+



MTA+ A Material for filling and rebuilding of root canals



Possible variants:

- inner dimension of the piston 0,8 mm inner dimension of the piston 1,2 mm

Reusable dental tools for handling and application of prepared measure of MTA+ material

MTA+ APPLICATOR A









Oil for dissolving gutta-percha

Orange oil 100%



Color chart and sizes:

ORANGE GUTTANE A





ENDO-ASPIRATOR A

for removal liquids from the root canal by aspirating









Choose the one according to your preferences

BLUE ETCH A

YELLOW

Dental etchant **0-phosphoric acid 36%** ▲ BLUE ETCH FLOW Flowable dental etchant 0-phosphoric acid 36%









Available packages: syringe 2 ml PORCELAIN ETCH A Porcelain etch Hydrofluoric acid 9,5%



▲ SILANE Methacrylic silane

24 RESTORATIVE DENTISTRY





SPECIAL



ALUSTAT GEL GEL for staunching the bleeding Aluminium chloride 25%



ALUSTAT LIQUID for staunching the bleeding Aluminium chloride 25%

Available package: bottle of 10 g





RESTORATIVE DENTISTRY 27



Available packages: dropping bottle 10 ml, 20 ml

▲ EUGENOL Eugenol oil Eugenol oil 100%







▲ ZINC OXIDE Powder for making pastes ZINC OXIDE Zinc oxide - fast setting formula



▲ RED DETECTOR Caries detector Rhodamine








PULP SPRAY SPRAY for pulp vitality testing















▲ RUBBER-DAM CORD Dental dam stabilizing cord



new

applicators






For:

liquids • creams • pastes • gels

Flexibility and right bending angle allows the proper application of preparations. The applicators have a universal luer-lock thread.



Elastic applicator 27G

lenght: **35** mm diameter: **0,4** mm



lenght: **41** mm diameter: **0,4 - 0,2** mm



lenght: **41** mm diameter: **0,4 -0,2** mm



Elastic applicator 25G

lenght: **41** mm diameter: **0,5 - 0,3** mm



lenght: **56** mm diameter: **0,5 -0,3** mm

Elastic applicator 22 G

lenght: **32,3** mm diameter: **0,7 - 0,5** mm

Elastic applicator 21 G

lenght: **56** mm diameter: **0,8 -0,6** mm







ADDITIONAL ACCESSORIES

Available packages:

1. ELASTIC APPLICATORS 10 pcs or 20 pcs

2. GEL APPLICATORS 20 pcs or 50 pcs

3. BRUSH APPLICATORS 10 pcs or 20 pcs

4. APPLICATORS FOR ALUSTAT FOAM 10 pcs

5. APPLICATORS FOR SYNTEX (mixing tips) 5 pcs

6. PAPER FOR PRODUCT MIXING 30 pcs

7. APPLICATION PIPETTES 5 pcs










ENDO-STATION GO

ENDO-STATION

Line of endodontic organizers









ENDO-STATION GO available in colors:*

green blue yellow orange orange neon transparent

* ENDO-STATION GO GUTTA BOX available only transparent

ENDO-STATION GO GUTTA BOX 2X3 ▶



ENDO-STATION GO MINI

ENDO-STATION GO MINI



ENDO-STATION GO MINI available in colors:

green blue transparent



- universal
- single use
- flexible



DENTAL CONTRASTOR





Perfect solution:

- V In a trip
- V For a business meeting
- V Before a date
- V During work
- V After every meal





Dent Fresh Smart 12 x 15 ml



Dent Fresh Smart KIT 10x15 ml + tooth cleaners small + tooth cleaners big





DENT FRESH

Two variants available:

- oryginal

- Available packages: Start Pack concentrate in syringe 50 ml
- + bottle
- + dispenser
- concentrate in syringe 50 ml









AVAILABLE PACKAGES:

Dent Fresh Smart Toothpaste 14 x 2 ml

DENT FRESH SMART

Toothpaste Original





Endo Removal System UPGRADED by K. Gończowski DDS, PhD, MSc

The ERS system is a universal set of instruments dedicated to the removal of broken parts of endodontic instruments from root canals under visual control by surgical microscope.

All parts of ERS are designed for multiple use



017-008

Endodontic

micro lever

www.endo-rs.com

015-195

Endodontic

micro probe - bent















PPH CERKAMED Wojciech Pawłowski

Address: ul. Kwiatkowskiego 1, 37-450 Stalowa Wolawoj. podkarpackie, Poland

Phone: +48 15 840 08 09

Phone: +48 15 842 35 85 in

Phone: +48 15 842 35 85, int. 809

Mobile: +48 601 278 010

E-mail: export@cerkamed.pl **WWW:** www.cerkamed.com

VAT: PL 8652048770 **REGON:** 431162810